



Clear revolution in orthodontics



Give Your Patients More Reasons To Smile

- ✓ CLEAR
- ✓ REMOVABLE
- ✓ COMFORTABLE
- ✓ HYGIENIC
- ✓ SNUG FIT
- ✓ DIGITAL DIAGNOSTIC

✉ support@dentments.com

🌐 www.dentments.com

The smile you deserve....



DentMents introduces an exciting new improved quality aligner through its unique process. Our mission is to deliver high quality, affordable, clear orthodontic aligner solution to the global market.

DentMents aligner provides a hygienic, convenient and clear solution to the correction of malocclusion without having to wear brackets and wires. Our staff members have over seven years of experience in the field of clear orthodontics. **DentMents** has completed a number of successful clinical trial cases before starting its commercial cases and has proven to be equally

successful in correction of a wider range of malocclusion.

DentMents Aligner and why is it is new?

- **DentMents** aligners are removable medical grade plastic appliances which patients wears instead of brackets and wires to correct malocclusion.
- A clear change to the concept of adult orthodontic treatment without using conventional methods of wires and brackets
- Adults wouldn't have asked for a better choice to correct those W4 malocclusions without having to

think about the social inhibitions.

- Predictable result visualization even before the treatment begins.
- Comfort to keep living a normal routine, without even modifying the eating habits.
- Hygienic choice for the correction of malocclusion, since it is clear and smooth all round, this makes it easier to maintain oral hygiene and its removable feature makes it the best hygienic choice.
- You are in control of treatment plan at every step of the way and can make suggestions for efficient flow of treatment to achieve desirable results.



DentMents Aligners



Metal wires and brackets



....is
within
your
reach
now

Why DentMents Aligners?

- No metal wires and brackets required to correct malocclusion
- Super fit, accuracy of movements
- Nearly invisible and clearer aligners
- See treatment outcome before even starting it
- No social inhibitions for adult patients to opt for orthodontic treatment
- Less chair side adjustment, shorter appointments
- More practice benefits

How to get DentMents Aligners? (case processing)

- Get registered with **DentMents**
- Gather patient records:-

Impression:

both upper and lower impression in rubber base impression material or STL Files

Photos:

Intra oral and extra oral profile photographs of patient

Radiograph:

Cephalometric and Panorex

Diagnostic and treatment form:

Filled completely


- Submit all the records to **DentMents** orthodontics
- Proposed diagnostic setup posted on the web based application for approval
- Aligner manufacturing: After approval of diagnostic setup aligners are manufactured
- Aligners are shipped to the doctor's office
- Patient wears aligners for a given period of time to achieve the desired smile.

BEFORE



AFTER





DentMents
CLEAR ALIGNERS

FOLLOW US AT



DENTMENTS



support@dentments.com



www.dentments.com